

**Wisconsin Oven Distributors, LLC.**

W355 S9075 Godfrey Ln.,

Eagly, WI 53119

Phone: 262.594.3941 • Fax: 262.594.3939

E-Mail: [tslaboven@centurytel.net](mailto:tslaboven@centurytel.net) • Web site: [www.memmertbenchlabovens.com](http://www.memmertbenchlabovens.com)

## Memmert CO<sub>2</sub> Incubators



CO<sub>2</sub> Incubator INCOmed  
"Celsius" Standard Software

Model Sizes: 108/153/246  
+20 °C to +50 °C  
Humidity 88 to 97 % rh  
CO<sub>2</sub> content up to 10 %  
Humidity 40 to 97 % rh (humidity module)  
CO<sub>2</sub> content up to 20 % (CO<sub>2</sub> module)  
O<sub>2</sub> content between 1 and 20 % (O<sub>2</sub> module)

**CO<sub>2</sub> INCUBATOR INCOmed** For cell cultivation and especially for in-vitro fertilisation, the precision and reliability of CO<sub>2</sub> incubators are of crucial importance. During cultivation, the slightest deviation in the CO<sub>2</sub> atmosphere, temperature or humidity can influence cell development. For this reason, Memmert has subjected its CO<sub>2</sub> incubators to a comprehensive evaluation process for their recognition as medical devices. The classification as class IIa medical device confirms that all Memmert CO<sub>2</sub> incubators INCOmed comply with the basic safety requirements of the European Medical Devices Directive 93/42/EEC. The interior chamber including all installations and sensors can be sterilised at 160 °C in a 4-hour programme.

Cultivation of Cells or Tissue, In-Vitro Fertilisation, Gene Expression

---

Results 1 - 3 of 3

Model	Previous Models	Approximate Volume	Width (A)	Height (B)	Depth (Less 25 mm for Fan) (C)
INCO108MED	INCO108	Approx. 3.8 ft <sup>3</sup> 108 L	22.0 in 560 mm	18.9 in 480 mm	15.7 in 400 mm
INCO153MED	INCO153	Approx. 5.4 ft <sup>3</sup> 153 L	18.9 in 480 mm	25.2 in 640 mm	19.7 in 500 mm
INCO246MED	INCO246	Approx. 8.7 ft <sup>3</sup> 246 L	25.2 in 640 mm	25.2 in 640 mm	23.6 in 600 mm

---

Results 1 - 3 of 3