

## CONSUMABLES

### for gas analyser OXYBABY® / PA / MAPY / MFA



Your WITT-Gas analyser is made of robust material and needs very little maintenance. There is a natural wear in the daily work. To make sure that replacement materials harmonize ideally with your WITT device we recommend original WITT spare parts. Cost-effective quick and reliable.

Suitable for all OXYBABY®, PA-, MAPY- and MFA- models.



#### O<sub>2</sub> sensor

Especially long life time (approx. 2 years in air). Can be replaced by customer self.



#### Filter



#### Needle



#### Foam rubber seals for sealing the sample site at the package

Package each 500 pieces



#### Hose with Luer-Lok connection

for using OXYBABY® as stationary device (table top version, with Can-Piercer or with Aquacheck)

